

Fabrikant apparaat	Type CPAP-apparaat	artikelnaam Filter	Filter type	Frequentie filter vervangen	Jaarvoorraad verstuurd in onderhoudsset
Breas	Isleep 20/20+/20i/22/25	Filter Stof, Isleep 20/20+/20i/22/25	HVF04-BR-FINE of 7270-407	1 x per 6 maanden	2 filters
Breas	PV-10	Filter Stof, PV-10	001972	1 x per 6 maanden	2 filters
Breas	Isleep 20/20+/20i/22/25	Filter Fijn Stof, Isleep 20/20+/20i/22/25	7270-399 of 004153 of	1 x per maand	12 filters
Covidien	Sandman info	Filter Stof, COV Sandman info	m41484006 (7270-450)	1 x per 6 maanden	2 filters
Covidien	Goodknight 420e/420s/425	Filter Stof, COV Goodknight 420e/420s/425	m41395002 (7270-252)	1 x per 6 maanden	2 filters
Fisher & Paykel	Icon series	Filter Stof, F&P Icon series	900icon503 (7270-499)	1 x per 3 maanden	4 filters
Fisher & Paykel	600 series	Filter Stof, F&P 600 series	HCFF02-FINE	1 x per 3 maanden	4 filters
Hoffrichter	Hoffrichter Autotrend Flexline/Trend 110/BIPAP Trend 300	Filter Stof, HOFF apparatuur	HCFHR02-FINE/4	1 x per 3 maanden	4 filters
Philips	CPAP Harmony	Filter Stof, Harmony	1029330	1x per 6 maanden	2 filters
Philips	Remstar CPAP Auto, C-flex, Plus, Pro, BIPAP Pro 2	Filter Stof, Remstar	HCFR07-FOAM	1 x per 6 maanden	2 filters
Philips	Remstar M-serie, Remstar Philips, System One	Filter Stof, Philips/ Remstar M-serie, System One	HCFR06-FOAM	1 x per 6 maanden	2 filters
Philips	Remstar BIPAP Pro	Filter Stof, Remstar BIPAP Pro	622220	1 x per 6 maanden	2 filters
Philips	Remstar CPAP Auto, C-flex, Plus, Pro, BIPAP Pro 2	Filter Fijn Stof, Remstar	HCFR06-FINE/6	1 x per 2 maanden	6 filters
Philips	Remstar M-serie	Filter Fijn Stof, Remstar M-serie	HCFR08-FINE/6	1 x per 2 maanden	6 filters
Resmed	S8 Elite/Elite 2/Escape 2/Autospirit 2	Filter Stof, RM S8	1850005	1 x per 6 maanden	2 filters
Resmed	Autoset CS2	Filter Stof, CS2	26920-1	1x per 2 maanden	6 filters
Resmed	Autoset T	Filter Stof, RM Autoset T	1850002	1 x per 6 maanden	2 filters
Resmed	S7 Lightweight/Autoset Spirit/VPAP 3	Filter Stof, RM S7 Lightweight/Autoset Spirit/VPAP 3	7270-410#6 of 30907	1 x per 2 maanden	6 filters
Resmed	S6 Elite/Plus	Filter Stof, RM CPAP S6 Elite/S6 plus	21936 (7270-44)	1 x per 2 maanden	6 filters
Resmed	S9 Autoset/Elite	Filter Stof, RM S9 Autoset/Elite	HCFS04-FINE/2 of 7270-471	1 x per 6 maanden	2 filters
Somnetics	CPAP Transcend	Filter Transcend	503003	1x per 6 maanden	2 filters
Weinmann	Somno - comfort2(e)/balance(e)/soft2(e)	Filter Stof, WM Somno - Comfort2(e)/Balance(e)/Soft2(e)	HCFW03-FOAM	1 x per 6 maanden	2 filters
Weinmann	Somno - soft/smart2/vent	Filter Stof, WM Somno-Soft/Smart2/Vent	WM15428	1 x per 6 maanden	2 filters
Weinmann	Somnotron-4/e	Filter Stof, WM Somnotron-4/e	WM15429	1 x per 6 maanden	2 filters